

OTech Sample Kit Data Sheet

Product description

Sample kit includes 4 samples of printed silver layers (electrodes) on different surfaces that showcase the variety of our inks. The samples included are:

1. Transparent layer on PET
2. Non-transparent reflective layer on PET
3. Non-transparent layer on fabric
4. Non-transparent layer on paper

The sample sizes are around 3.5 cm x 4.5 cm or similar.



OTech sample kit box

Benefits

You can take out all samples to see how they look or to touch how they feel. You can use them for first simple tests like conductivity, flexibility, transparency (only on OTech T), and haze or more advanced tests like microscopy or surface analysis. Additionally, the casing of the sample kit includes details about our inks and the most important technical data for each sample, so you and your team can brainstorm your next innovation.

Directions for use and storage

You can use the samples however you want. However, it should be noted that it is very easy to leave fingerprints on the transparent slides. Tweezers or gloves are, therefore, advised. Store in a dry place (office conditions).



The inside of the sample kit box